



Fraction Collector C-660

Collect each drop of your precious sample

Reliably collecting the separated compounds is the ultimate goal of a preparative chromatography application. The Fraction Collector C-660 collects all your samples according to time, volume or peak mode in any type of tubes or vessels.



Reliable

Retrieve each fraction in the format you require



Smart

Precise peak detection for two detector signals



Adaptable

Free configuration of the collection area

C-660: Your most important benefits



Reliable

- Retrieve any fraction in the requested format
- Proven peak detection algorithm
- Large collection capacity for long runs or high flow rates
- Precise compound collection with programmable delay time
- Safe run with optional waste level sensor



Smart

- Choose your collection mode according to time, volume or peak
- Use up to eight detector signals
- Maximized purity of collected fraction with separation of peak front and tail
- Manual or computer control



Adaptable

- Teach-in function for customer specific racks (regular or irregular)
- Adjustable height for optimal integration in your workflow
- Choice of four collection patterns
- Multi language user interface

Complete your portfolio



Pump Module
C-601 / C-605
Pump



Sepacore® Flash
cartridges
Prepacked columns



UV-Vis Detector
C-640
UV-Vis detection



SepacoreControl
Software
Full control software

