

Model	ERD-03
No. of rows	6
Vessel capacity	300 ml
Reciprocating height	100 mm
No. of dip rate	5 to 50 dpm
Accuracy (at 37°C)	0.2
Stroke height	100 mm
Readable temperature range	15 to 50°C
Controlled temperature range	30 to 40°C
Temperature resolution	0.2°C
Display	Touch screen
Motor	Single shaft motor 10 kg cm
Heater	1000 Watts, 230 V AC, S.S. 304
Water level sensor	Float switch type level sensor
Material of construction	S.S. 316
Water bath	Clear acrylic moulded water bath
Bath capacity	44 litres

<b>Model</b>	<b>ERD-07</b>
No. of rows	6
Vessel capacity	300 ml
Reciprocating height	100 mm
No. of dip rate	5 to 50 dpm
Accuracy (at 37°C)	0.2
Stroke height	100 mm
Readable temperature range	15 to 50°C
Controlled temperature range	30 to 40°C
Temperature resolution	0.2°C
Display	Touch screen
Motor	Single shaft motor 10 kg cm
Heater	1000 Watts, 230 V AC, S.S. 304
Water level sensor	Float switch type level sensor
Material of construction	S.S. 316
Water bath	Clear acrylic moulded water bath
Bath capacity	44 litres