

Part No.	Volume Range (µl)	Volume Increments (µl)	Mode (M/RD) ¹	Test Volume (µl)	Accuracy ² (±%)	Precision ³ (≤%)
2011	0.2 - 2	0.001	M	0.2	12.0	6.0
				1	3.0	1.6
				2	1.5	0.7
			RD	0.5	9.0	9.0
2012	1 - 10	0.01	M	1	5.0	1.5
				5	1.5	0.8
				10	1.2	0.5
			RD	1	6.5	6.5
2013	2 - 20	0.01	M	2	7.0	2.0
				10	1.6	0.5
				20	1.0	0.4
			RD	2	4.5	4.5
2015	10 - 100	0.1	M	10	3.5	1.0
				50	1.0	0.35
				100	0.8	0.2
			RD	10	3.5	3.5
2016	20 - 200	0.1	M	20	2.5	1.0
				100	1.0	0.3
				200	0.7	0.2
			RD	20	3.0	3.0
2018	100 - 1000	1	M	100	3.0	0.6
				500	1.0	0.3
				1000	0.8	0.2
			RD	100	2.0	2.0

1. M: Manual mode, RD: Repeat Dispense mode

2. Accuracy: systematic error

3. Precision: random error, coefficient of variation